

## OTC




DriWeave® Woven Top Cover




NST® Ultra Thin Non-Skid Bottom



12mm Heel Cup Intrinsically Posted



3mm Polypropylene Shell With Intrinsic Forefoot Post



The team at Performance has developed a functional OTC orthotic that is only available through practicioner offices.

Linear Foot Motion devices use trusted podiatric designs to control your patients foot through out the gait cycle.

- Employing a 12mm heel cup and an intrinsically posted rear foot and forefoot 3mm polypropylene shell.
- A full foot EVA interface with woven fabric top cover to protect your patient from the forces of pressure, shock and sheer.
- An NST® non-slip extension bottom resists separation at the metatarsals and provides extended life to this beautifully crafted appliance that's the perfect addition to your practice.

\*\* Performance\*\*\*

Call Performance for further details.